

PRPP Assay Kit


**For a one-step enzymatic measurement of
 α -D-5-phosphoribosyl-1-pyrophosphate (PRPP)**

NOVO CIB's PRPP Assay Kit is designed to measure PRPP (α -D-5-phosphoribosyl-1-pyrophosphate) content in samples. This enzymatic assay is based on a coupled reaction involving Hypoxanthine-guanine phosphoribosyltransferase (HGPRT) and Inosine Monophosphate Dehydrogenase (IMPDH). PRPP is highly unstable, which renders the preparation of accurately weighted standards difficult. **NOVO CIB's** PRPP Assay Kit has the advantage of measuring PRPP concentration in the sample by absorbance through a stoichiometric NADH₂ formation, thus allowing a direct and absolute measurement of PRPP content with no need of any calibration curve.


Principle

NOVO CIB's PRPP Assay Kit involves the coupling of the following enzymatic reactions:

- (1) In the presence of Hypoxanthine, the 5-phosphoribosyl group of PRPP is first transferred by Hypoxanthine-guanine phosphoribosyltransferase (HGPRT) to form Inosine Monophosphate (IMP) and release a pyrophosphate:



- (2) IMP is immediately oxidized to Xanthosine Monophosphate (XMP) by a highly active IMPDH in the presence of NAD.



The NADH₂ formed is equivalent to the amount of PRPP present in the assay. NADH₂ formation is directly monitored spectrophotometrically at 340 nm.

Assay with a range of PRPP standards

Figure 1 shows the correlation between PRPP concentrations of standards ranging from 2 to 400 μ M in 200 μ l-well assays and the PRPP concentrations as measured spectrophotometrically at 340nm by applying a 0.625cm light-path (corresponding to 200 μ l/well in a standard 96-well microplate).

Kit Content

Reaction buffer x10

Reagents: NAD, Hypoxanthine

IMPDH enzyme, lyophilized *

HGPRT enzyme, lyophilized **

* Recombinant IMPDH from *Staphylococcus aureus*, expressed in *E. coli* (Ref. #E-Nov7)

** Human recombinant HGPRT, expressed in *E. coli* (Ref. #E-Nov9)




Figure 1: Correlation between PRPP concentration of the standards and PRPP concentration measured with NOVO CIB's PRPP Assay Kit.

Assays were done in duplicate using a range of freshly prepared PRPP (Sigma-Aldrich, P8296). The concentration of the standards indicated on the graphic takes into account the 75% purity indicated by the supplier.

Due to PRPP instability, standards may be difficult to prepare. This correlation curve shows that a direct measurement of PRPP content can be done without any standard but directly through the absorbance of the NADH₂ formed in the assay after 35 minutes of reaction.