

Development of *in vitro* and cell-based assays for assessing nucleotide biosynthesis inhibition

Larissa Balakireva, Blandine Gri, Thomas Liechti, Nicolas Godard NovoCIB, 115 avenue Lacassagne 69003 Lvon, France

www.novocib.com

Background Numerous nucleoside analogues (NAs) are currently used to treat viral hepatitis. They are usually designed to inhibit one viral target. This remains in contrast with the observation that ribavirin (1-β-D-ribofuranosyl-1,2,4-triazole-3-carboxamide), a purine nucleoside analogue currently used as a part of bi-therapy of hepatitis C infection, has multiple modes of action: (i) depletion of intracellular GTP pools by inhibition of the cellular IMPDH, (ii) inhibition of viral polymerase activity by the 5'-triphosphate metabolite of ribavirin, (iii) induction of error catastrophe as a result of accumulation of mutations in the viral genome. Even if no direct relationship between ribavirin antiviral action and IMPDH inhibition has been demonstrated, the depletion of cellular GTP should result in an increased frequency of ribavirin triphosphate incorporation by viral polymerase due to lower intracellular concentration of its natural competitor.

Aims Relation between therapeutic potential of a nucleoside analogue and its anti-metabolite action remains difficult to demonstrate mainly because of the lack of investigation tools. The purpose of this study was to develop a range of assays to reveal nucleotide biosynthesis inhibition.

Methods and Results For a rapid evaluation of new nucleoside analogues as IMPDH inhibitors, we have developed an in vitro enzymatic assay where the synthesis of monophosphorylated form of nucleoside analogue (NA-MP) is provided by cloned human nucleoside kinases, and NA-MP is immediately tested for inhibition of human recombinant IMPDH. This assay has been validated with nucleoside analogues ribavirin and mizoribine.

We have also developed original cell-based analytical approach in which 27 cellular ribo- and deoxyribonucleotides are extracted from cultured cells. separated by ion-pairing chromatography and quantified. This cellular assay, validated with several NA (ribavirin, aracytidine, gemcitabine) and known antimetabolites (mycophenolic acid, leflunomide, hydroxyurea), provides a powerful tool for studying the effect of new nucleoside analogues on whole spectra of cellular purine and pyrimidine deoxy- and ribonucleotides. In addition, in regards with new antiviral molecules identified in HCV cell culture systems, our cellular assay allows to distinguish the molecules that directly acts on the viral proteins from others that inhibits the cell nucleotide biosynthesis.




Ribavirin (1-b-D-ribofuranosyl-1,2,4-triazole-3-carboxamide, ref. 1) is a purine nucleoside analogue with broad spectrum antiviral activity. Since the 1970's, it is known that the initial step of ribayirin phosphorylation is provided by adenosine kinase (2). Recently it has been demonstrated that cytosolic 5'-nucleotidase II can also phosphorylate ribavirin, that could contribute to the development of ribavirin-induced haemolytic anemia in vivo (3). Monophosphorylated form of Rbv was shown to inhibit cellular IMPDH, whereas triphosphorylated form, after being incorporated by viral RNA-dependent RNA-polymerases, duces "error catastrophe" and inhibits RNA synthesis.

Enzymatic assays for the characterization of phosphorylation properties of new nucleoside analogues


The therapeutic efficacy of nucleoside analogues is highly dependent on their intracellular phosphorylation. Nucleoside phosphorylation assays focused on cloned human recombinant nucleoside kinases (dCK, AK, YMPK, cNII) have been developed and validated with known NAs (Table 1). Vmax and Km values of human recombinant nucleoside kinase were determined for natural substrates and nucleoside analogues (e.g. ribavirine, araC, gemcitabine, etc). The results obtained are highly similar to previously published data. These assays now allow the characterization of substrate properties (Km and Vmax) of new nucleoside analogues in comparison with those of known nucleoside analogues

Table 1. Cloned human recombinant nucleoside kinases and their substrates.



Enzymatic characterization of human recombinant IMPDH II

IMP Dehydrogenase (IMPDH, E.C. 1.1.1.205) catalyzes the pivotal step in guanine nucleotide biosynthesis: the conversion of IMP to XMP. Blocking the conversion of IMP to XMP. Rbv leads to depletion of GTP pools, that results in increased frequency incorporation of Rbv-TP by viral polymerase due to intracellular concentration of its natural



Enzymatic assay with cloned human IMPDH II was developed and validated by measuring Vmax and Km value for IMP and IC50 value for MPA.


Cellular nucleoside kinases

 $IMP \longrightarrow XMP \longrightarrow GMP \longrightarrow GDP \longrightarrow GTP$


viral

- [1] D.G. Streeter, J. T. Witkowski, G.P. Khare, R.W. Sidwell, r.j. Bauer, R.K. Robins, L. N. Simon (1973) Mechanism of Action of 1-t3-D-Ribofuranosyl-1,24-friazole-3-Carboxamide (Virazole), A New Broad-Spectrum Antiviral Agent Proc. Nat. Acad. Sci. USA Vol. 70, No. 4, pp. 1174-1178
- R. C. Wills, D. A. Carson, and J. E. Seegmiller (1973) Adenosine Kinase Initiates the Major Route of Ribavirin Activation in a Cultured Human Cell Line. PNAS USA, 75: 3042-3044.
 Wu JZ, Larson G, Walker H, Shim JH, Hong Z. Phosphorylation of ribavirin and viramidine by adenosine kinase and cytosolic 5-nucleotidase II: Implications for ribavirin metabolism in erythrocytes. (2005) Antimicrob Agents Chemother. 49(6):2164-71

Coupled nucleoside kinase -IMPDH II assay for rapid evaluation of monophosphate forms of nucleoside analogues




phosphorylated to inhibit IMPDH. For rapid evaluation of NA-monophosphate as IMPDH enzymatic assay based on cloned human as ADP formed B. The IMPDH action



In this assay, the monophosphorylation step of nucleoside analogue is provided by specific nucleoside kinases; the enzymatically produced NA-monophoshate is then directly tested for IMPDH inhibition without puriifcation. This combined test has been validated with Rby and mizoribine.

Whole cell bio-analytic assays: nucleotide profiling

To allow the study of whole spectra of purine and pyrimidine ribo- and deoxyribonucleotides in drug-treated cultured cells, we have developed whole cell bio-analytic assay. More than 27 (deoxy)ribonucleotides (mono-, di-, triphosphate) and nucleotide co-factors are extracted from cultured cells, separated by ion-pared chromatography and quantified. This cellular assay was validated both with anti-viral and anti-cancer nucleoside analogues and known



In conclusion, we have developed a range of new tools that allow

- 1) to characterize in vitro phosphorylation properties of new nucleoside analogues, using recombinant human nucleoside kinases
- 2) to characterize IMPDH inhibition by monophosphorylated form of NAs
- 2) to study anti-metabolite effect of new antiviral molecules and to distinguish "true" antiviral molecules from inhibitors of cell

